

# NObreath<sup>®</sup>

## Mouthpiece Datasheet



Designed specifically for the NObreath<sup>®</sup> FeNO monitor, the NObreath<sup>®</sup> mouthpiece contains an integrated infection control filter that removes and traps <99% of airborne bacteria and >98% of viruses<sup>1</sup>.

### Technical Specification

Dimensions	Approx. 180 x 28 x 22 mm
Weight	Approx. 14g
Material	Polypropylene
Shelf life	5 Years



The NObreath<sup>®</sup> mouthpiece may be used multiple\* times by one patient, per clinical visit.

**NOTE:** The mouthpiece is single-patient use; reuse of the mouthpiece between patients can increase the risk of cross-infection and result in incorrect readings.

The NObreath<sup>®</sup> mouthpiece is fully recyclable and should be disposed of after use in accordance with local clinical waste guidelines.

#### References:

1. Public Health England. An Evaluation of Filtration Efficiencies Against Bacterial and Viral Aerosol Challenges. Salisbury: Public Health England; 2020

\*The mouthpieces are single patient use only and can be used for a maximum of 3 tests

Visit [www.bedfont.com/resources](http://www.bedfont.com/resources) to view this document in other languages.



Bedfont<sup>®</sup> Scientific Ltd.  
Station Road, Harrietsham, Maidstone,  
Kent, ME17 1JA England  
Tel: +44 (0)1622 851122 Fax: +44 (0)1622 854860  
Email: [ask@bedfont.com](mailto:ask@bedfont.com) Web: [www.bedfont.com](http://www.bedfont.com)



Stephen Rowe  
Cristimar E4-1  
Ave Juan Carlos I  
Los Cristianos, Arona, 38650  
Santa Cruz de Tenerife, Spain

© Bedfont<sup>®</sup> Scientific Limited 2022

Issue 6 - May 2022, Part No: MKT514

Bedfont<sup>®</sup> Scientific Limited reserves the right to change or update his literature without prior notice. Registered in: England and Wales. Registered No: 1289798

