

NObreath[®]

Mouthpiece Datasheet



Designed specifically for the NObreath[®] FeNO device, the NObreath[®] mouthpiece contains an integrated infection control filter that removes and traps >99% of airborne bacteria and >98% of viruses¹.

Technical Specification

Dimensions	Approx. 180 x 28 x 22 mm
Weight	Approx. 14g
Material	Polypropylene
Shelf life	5 Years

WARNING: The mouthpieces are single patient use only and can be used for a maximum of 10 tests per breath testing session. Further re-use could cause incorrect readings and could increase the risk of cross infection. The mouthpiece should be disposed of after use, in accordance with local waste disposal guidance.



References:

1. Public Health England. An Evaluation of Filtration Efficiencies Against Bacterial and Viral Aerosol Challenges. Salisbury: Public Health England; 2020.

Visit www.bedfont.com/resources to view this document in other languages.



Bedfont[®] Scientific Ltd.
Station Road, Harrietsham, Maidstone,
Kent, ME17 1JA England
Tel: +44 (0)1622 851122 Fax: +44 (0)1622 854860
Email: ask@bedfont.com Web: www.bedfont.com



Emergo Europe B.V.
Westervoortsewijk 60
6827 AT Arnhem
The Netherlands.

© Bedfont[®] Scientific Limited 2024

Issue 9 - September 2024, Part No: MKT514

Bedfont[®] Scientific Limited reserves the right to change or update this literature without prior notice. Registered in: England and Wales. Registered No: 1289798



MD 502805